


Exceed Your Expectations





HEMODIALYZER
REXEED™
REXBRANE™

Advanced Engineering




Dye was injected into the dialysate inlet port, with dialysate flow at 500 mL/min, then 10mL fractions were collected at the outlet port and measured for dye concentration by absorbance.



A sharp, symmetrical peak indicates rapid and complete penetration of dye into the center of the fiber bundle to give uniform flow distribution, proving the benefits of fiber bundle and housing design optimization.

Therapeutic Benefits

Clinical Evaluation of REXEED-A and Other PS




Good phosphate management, particularly in higher body weight patients is afforded by the superior clearance characteristics of REXEED-A.

Maurizi J. et al., CHU (University Hospital Center) Michallon, Grenoble, France. In vivo Clinical Evaluation, 2005.


REXEED-A has high removal for small toxins as well as excellent reduction rates for low molecular weight proteins such as β_2 -MG (M.W. 11,800). The sharp membrane cut-off minimizes albumin loss.

Clinical Evaluation of REXEED-A and PES



Adapted from: S. Osawa, H. Yamamoto, S. Hisajima, "Evaluation of dialyzers with various characteristics", Jpn. J. High Performance Membrane '05, Kidney & Dial. Suppl. 59: 216-219, 2005.

Pre - Dialysis Creatinine



Three high body mass patients were changed from their original high flux polysulfone dialyzers (2.2m²) to REXEED-25A and dialyzed for 4 weeks. Dialysis conditions were unchanged during the evaluation period, after which the patients were returned to their original dialyzers.

Lindley E.J. et al., Leeds Teaching Hospitals NHS Trust, Leeds, UK. In vivo Clinical Evaluation, 2005.

REXEED

The high performance dialyzer for high performance dialysis

REXEED™ - A
Performance (in vitro)

		REXEED-13A [※]			REXEED-15A			REXEED-18A			REXEED-21A			REXEED-25A		
Blood Flow (mL/min)		200	300	400	200	300	400	200	300	400	200	300	400	200	300	400
Clearances (mL/min)																
	Urea	194	264	306	196	272	317	198	280	330	199	284	340	199	287	350
	Creatinine	188	247	284	191	255	296	194	265	309	195	272	317	196	280	330
	Phosphate	180	227	259	185	238	272	190	250	289	192	257	299	194	260	303
	Vitamin B ₁₂	133	155	165	142	168	179	152	183	197	160	194	212	163	199	225
	Inulin (Qb=200mL/min)	78			87			92			102			117		
KoA (Urea) (mL/min)		1027			1190			1415			1569			1714		
KUF (mL/hr/mmHg (mL/hr/kPa))		66 (495)			72 (540)			81 (608)			90 (675)			102 (773)		
Effective Surface Area (m ²)		1.3			1.5			1.8			2.1			2.5		
Sieving Coefficient																
	Inulin	1.0														
	β ₂ -microglobulin	0.8														
	Myoglobin	0.4														
	Albumin	0.001														

Conditions: Clearances: Qd=500mL/min, Qf=0mL/min KoA: Qb=300mL/min, Qd=500mL/min, Qf=0mL/min
KUF: Measured with bovine blood, TP=60±5g/L, Hct=32±3%, Qb=300mL/min
ISO8637-1:2017

Specifications

	REXEED-13A [※]	REXEED-15A	REXEED-18A	REXEED-21A	REXEED-25A
Membrane	REXBRANE™ (Asahi Polysulfone)				
Internal Diameter of Hollow Fiber (μm)	185				
Wall Thickness of Hollow Fiber (μm)	45				
Priming Volume (mL)	76	86	103	117	137
Maximum TMP (mmHg (kPa))	500 (66)				
Maximum Blood Flow (mL/min)	500				
Maximum Dialysate Flow (mL/min)	800				
Dimensions (mm [L] x mm [D])	334 x 38	334 x 41	334 x 43	334 x 47	334 x 50
Weight (g)	412	463	520	604	675
Sterilization	Gamma-Ray				

Note: High permeability devices. Use only with ultrafiltration controlling equipment. ISO8637-1:2017
* Asia limited countries only

AsahiKASEI

Manufacturer:
ASAHI KASEI MEDICAL CO., LTD.
1-1-2 Yurakucho, Chiyoda-ku, Tokyo 100-0006 Japan
TEL: +81-3-6699-3771 FAX: +81-3-6699-3773
www.asahikasei-medical.com

Represented in Europe by:
ASAHI KASEI MEDICAL EUROPE GmbH
Herriotstrasse 1, 60528 Frankfurt am Main Germany
TEL: +49 (0) 69 66 37 15 00 FAX: +49 (0) 69 6 66 51 93

