
# LIPASE [LP]

from *Chromobacterium viscosum*  
 (Triacylglycerol acylhydrolase, EC 3.1.1.3)  
 (Triacylglycerol lipase)

## Preparation and Specification


Appearance	: White to off-white amorphous powder, lyophilized
Specific activity	: More than 2,500 U/mg solid
Contaminants	:
Cholesterol oxidase	Less than 0.01 % (U/U)
Catalase	Less than 0.01 % (U/U)

## Properties

Substrate specificity	: See Table 1	
Molecular weight	: 120 kDa (pH 3.7) (gel filtration) 30 kDa (pH 7.3) (gel filtration)	
Isoelectric point	: pH 3.7 and pH 7.3	
Optimum pH	: 3.0–10.0	Figure 1
pH stability	: 4.0–10.0 (50°C, 60 min)	Figure 2
Thermal stability	: Stable at 70°C and below (pH 7.0, 10 hr)	Figure 3
Storage stability	: At least one year at -20°C	Figure 4
Effect of various chemicals	: See Table 2	
Activators	: Detergents	

## Applications for Diagnostic Test

This enzyme is useful for enzymatic determination of **triglyceride** when coupled with glycerophosphate oxidase (T-60) and glycerol kinase (T-64).

TG: Triglyceride

FFA: Free fatty acid

DHAP: Dihydroxyacetone phosphate


Table 1. Substrate specificity

Substrate	Relative activity (%)
Triolein	100
Tripalmitin	22
Trimyristin	53
Trilaurin	103
Tricaprin	166
Tricaprylin	312
Tricaproin	156
Tributylin	94
Tripropionin	22
Triacetin	38

Table 2. Effect of various chemicals on LP activity


Additives	Concentration	Relative activity (%)
None	-	100
NiCl	1mM	97
MnCl <sub>2</sub>	1mM	98
(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	1mM	100
MgCl <sub>2</sub>	1mM	99
ZnCl	1mM	94
ZnSO <sub>4</sub>	1mM	94
Ba(CH <sub>3</sub> COO) <sub>2</sub>	1mM	98
CaCl <sub>2</sub>	1mM	100
MoSO <sub>4</sub>	1mM	106
CuSO <sub>4</sub>	0.5mM	26
CuCl <sub>2</sub>	0.5mM	26
FeCl <sub>3</sub>	1mM	103
CoCl <sub>2</sub>	1mM	96
Li <sub>2</sub> CO <sub>3</sub>	1mM	93
EDTA	1mM	105
KCl	100mM	100
NaCl	100mM	98
NaN <sub>3</sub>	0.05%	97
NaF	20mM	88

Fig.1 pH Optimum


▲ : Mcllvaine buffer  
 □ : Phosphate buffer  
 △ : Borate buffer

Fig.2 pH Stability


● : 37°C, 24 hr  
 ○ : 50°C, 1 hr  
 pH 3, 4, 5 Mcllvaine buffer  
 pH 6, 7, 8 Phosphate buffer  
 pH 9, 10, 11 Borate buffer

Fig.3 Thermal Stability

○ : Lyophilized powder  
 ● : Aqueous solution

Fig.4 Storage (lyophilized powder)


○ : -20°C  
 □ : 5°C  
 △ : 30°C


## Assay

### Principle

The assay is based on the titration of fatty acids liberated in the following reactions:

(Titration)

### Unit definition

One unit is defined as the amount of enzyme which liberates 1 μ mole of fatty acid per minute at 37°C under the conditions specified in the assay procedure.

### Reagents

- Substrate suspension (Olive oil and Adekatol SO-120)  
 50 g of Olive oil (Japanese Pharmacopoeia grade) and 50 g of Adekatol SO-120 are suspended with 150 ml of distilled water.

### 2. Reaction stopper

Mixture of ethanol and acetone (1:1)

### 3. Indicator

1% (W/V) Phenolphthalein-ethanol solution

### 4. Titration solution

50 mM NaOH solution

### 5. Enzyme dilution buffer

0.1 M KH<sub>2</sub>PO<sub>4</sub>-NaOH buffer, pH 8.0 containing 0.1% (W/V) BSA and 0.1% (W/V) NaN<sub>3</sub>

### 6. Reagents

Olive oil: (Japanese Pharmacopoeia grade)

Ethanol: (Japanese Pharmacopoeia grade)

Adekatol SO-120 : ADEKA CORPORATION

BSA: Millipore Fraction V pH5.2 #81-053

### Enzyme solution

Accurately weigh about 10 mg of the sample and add enzyme dilution buffer to make a total of 50 ml.

Dilute it with enzyme dilution buffer to adjust the concentration to within 2-4 U/ml.

## ■ Procedure

1. Pipette accurately 5 ml of substrate suspension and 2 ml of distilled water into a test tube (24 mm i.d. × 200 mmH) and mix to start the preincubation at 37°C.
2. After 10 min, add 0.5 ml of enzyme solution and mix to start the reaction.

※ In the case of a test blank, add 0.5 ml of enzyme dilution buffer in place of enzyme solution.

3. After 20 min, stop the reaction with 16 ml of reaction stopper.
4. Add 3 drops of indicator and titrate the whole mixture with under nitrogen gas bubbling.

※ End point of titration: Appearance of the same color as that of the blank

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20 : reaction time (min)

F : factor of titration solution (50 mM NaOH)

50 : concentration (mM) of titration solution (50 mM NaOH)

0.5 : the volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

## Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition (Figure 4).

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