
(Diagnostic Reagent Grade)

ASAHI KASEI ENZYMES T-222

PHOSPHOLIPASE D [PLD II]

from *Streptomyces chromofuscus*
(Phosphatidylcholine phosphatidohydrolase: EC 3.1.4.4)

Preparation and Specification

Appearance : Pale grayish to grayish or brownish to light purple lyophilizate
 Specific activity : More than 30 U/mg solid
 Contaminants :
 Catalase : Less than 0.6 % (U/U)
 Glucose oxidase : Less than 0.02 % (U/U)

Properties

Substrate specificity : See Table 1
 Molecular weight : 58 kDa (SDS-PAGE)
 Isoelectric point : pH 5.9 (estimated from amino acid sequence)
 Michaelis constants : 1,2-Dioleoyl-sn-glycero-3-phosphocholine 9.3×10^{-4} M
 Optimum pH : 6.7-7.1 Figure 1
 pH stability : 5.3-9.7 Figure 2
 Thermal stability : Stable at 60°C and below (pH 8.0, 10 min) Figure 3
 Storage stability : At least one year at -20°C
 Effect of metal ions : See Table 2
 Activators : Ca²⁺

Applications for Diagnostic Test

This enzyme is useful for enzymatic determination of **phospholipids** when coupled with choline oxidase (T-05)



Table 1. Substrate specificity

Substrate	Specific activity (%)
1,2-Dioleoyl-sn-glycero-3-phosphocholine	100
2-Oleoyl-1-palmitoyl-sn-glycero-3-phosphocholine	98
L- α -Phosphatidylcholine	95
L- α -Lysophosphatidylcholine	99
1-Oleoyl-sn-glycero-3-phosphocholine	99
L- α -phosphatidylethanolamine	14
Sphingomyelin	26

Table 2. Effect of metal ions (Activators)


Divalent cation (1 mM)	Relative activity (%)
None	2
Ca ²⁺	100
Mg ²⁺	2
Mn ²⁺	0
Ba ²⁺	1

Fig.1 pH Optimum


○ : 3,3-Dimethylglutarate-NaOH buffer
 □ : Tris-HCl buffer
 ● : Glycine-NaOH buffer

Fig.2 pH Stability

37°C, 60 min.
 10 mM buffer containing 0.1% TritonX-100 and 0.05% BSA
 ○ : 3,3-Dimethylglutarate-NaOH buffer
 □ : Tris-HCl buffer
 ● : Glycine-NaOH buffer

Fig.3 Thermal Stability


pH 8.0, 10 min.
 10 mM Tris-HCl buffer containing 0.1% TritonX-100 and 0.05% BSA

Assay

Principle

The assay is based on the increase in absorbance at 500 nm as the formation of quinoneimine dye proceeds in the following reactions:

COD: Choline oxidase

Unit definition

One unit is defined as the amount of enzyme which hydrolyzes 1 μ mole of phosphatidylcholine to phosphatidic acid and choline per minute at 37°C under the conditions specified in the assay procedure.

Reagents

- Reaction mixture for the first reaction

0.1 M Tris-HCl buffer pH 8.0	0.20 ml
0.1 M CaCl ₂ solution	0.05 ml
25 mM substrate solution ¹⁾	0.10 ml

- | | |
|------------------------------|---------|
| Distilled water | 0.15 ml |
| 1): 25 mM substrate solution | |

Dissolve 88.5 mg of 1,2-Dioleoyl-sn-glycero-3-phosphocholine with 4.5 ml of 5 % (W/V) Triton X-100 solution.

- Reaction mixture for the second reaction

15 mM 4-AA solution	0.10 ml
0.2 % (W/V) Phenol solution	0.10 ml
60 mM EDTA pH 8.0	0.10 ml
50 mM Tris-HCl buffer pH 8.0	2.00 ml
90 U/ml POD solution ²⁾	0.10 ml
30 U/ml COD solution ³⁾	0.10 ml

EDTA: Ethylenediaminetetraacetic acid

2): 90 U/ml POD solution

Dissolve 900 U (PPU) of POD with 10 ml of distilled water.

3): 30 U/ml COD solution

Dissolve 300 U of COD with 10 ml of 10 mM Tris-HCl buffer pH 8.0.

- Enzyme dilution buffer

10 mM Tris-HCl buffer (pH 8.0) containing 0.05% (W/V) BSA and 0.1% (W/V) Triton X-100

4. Reagents

Triton X-100: The Dow Chemical Company

1,2-Dioleoyl-sn-glycero-3-phosphocholine:

Sigma Chemical Co. #P-6354

EDTA (2 Na·2H₂O): KISHIDA CHEMICAL Co., Ltd.

#060-29133

COD: Asahi Kasei Pharma Corporation #T-05
 BSA: Millipore Fraction V pH 5.2 #81-053
 4-AA: NACALAI TESQUE, INC. Special grade #01907-52
 POD: Sigma Chemical Co. Type II #P-8250

Absorbance sample : As
 blank : Ab
 $\Delta A = (As - Ab) \leq 0.60 \text{ Abs}$

Enzyme solution

Accurately weigh about 20 mg of the sample and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

Procedure

- Pipette accurately 0.50 ml of reaction mixture for the first reaction into a small test tube and preincubate at 37°C.
- After 5 min, add 50 μ l of enzyme solution and mix to start the reaction at 37°C.
- At 10 min after starting the reaction, add 2.50 ml of reaction mixture to the second reaction and mix to start the second reaction.
 ※ In the case of a blank test, add 50 μ l of enzyme dilution buffer solution at this time.
- At 20 min after starting the reaction, measure the absorbance at 500 nm.

Calculation

$$\text{Activity (U/mg of powder)} = \frac{\Delta A / 10}{12.2} \times \frac{3.05}{0.05} \times \frac{1}{X}$$

12.2 : millimolar extinction coefficient of quinoneimine dye
 ($\text{cm}^2 / \mu\text{mole}$)

10 : reaction time (min)

3.05 : final volume (ml)

0.05 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution
 (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition.

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