
Table 1. Substrate specificity


Substrate	Relative activity (%)
Triolein	100
Tripalmitin	22
Trimyristin	53
Trilaurin	103
Tricaprin	166
Tricaprylin	312
Tricaproin	156
Tributylin	94
Tripropionin	22
Triacetin	38

Fig.1 pH Optimum


▲ : Mcllvaine buffer
□ : Phosphate buffer
△ : Borate buffer

Fig.2 pH Stability


● : 37°C, 24 hr
○ : 50°C, 1 hr
pH 3, 4, 5 Mcllvaine buffer
pH 6, 7, 8 Phosphate buffer
pH 9, 10, 11 Borate buffer

Fig.3 Thermal Stability

pH 7.0, 10 hr
○ : Lyophilized powder
● : Aqueous solution

Fig.4 Storage (lyophilized powder)


○ : -20°C
□ : 5°C
△ : 30°C

Assay

Principle

The assay is based on the titration of fatty acids liberated in the following reactions:

Unit definition

One unit is defined as the amount of enzyme which liberates 1 μ mole of fatty acid per minute at 37°C under the conditions specified in the assay procedure.

Reagents

- Substrate suspension (Olive oil and Adekatol SO-120)
50 g of Olive oil (Japanese Pharmacopoeia grade) and 50 g of Adekatol SO-120 are suspended with 150 ml of distilled water.
- Reaction stopper
Mixture of ethanol and acetone (1:1)
- Indicator
1% (W/V) Phenolphthalein-ethanol solution

4. Titration solution

50 mM NaOH solution

5. Enzyme dilution buffer

0.1 M KH_2PO_4 -NaOH buffer, pH 8.0 containing 0.1% (W/V) BSA and 0.1% (W/V) NaN_3

6. Reagents

Olive oil: (Japanese Pharmacopoeia grade)
Ethanol: (Japanese Pharmacopoeia grade)
Adekatol SO-120: ADEKA CORPORATION
BSA: Millipore Fraction V pH5.2 #81-053

Enzyme solution

Accurately weigh about 10 mg of the sample and add enzyme dilution buffer to make a total of 50 ml.

Dilute it with enzyme dilution buffer to adjust the concentration to within 2-4 U/ml.

Procedure

- Pipette accurately 5 ml of substrate suspension and 2 ml of distilled water into a test tube (24 mm i.d. \times 200 mmH) and mix to start the preincubation at 37°C.
- After 10 min, add 0.5 ml of enzyme solution and mix to start the reaction.
* In the case of a test blank, add 0.5 ml of enzyme dilution buffer in place of enzyme solution.

- After 20 min, stop the reaction with 16 ml of reaction stopper.
- Add 3 drops of indicator and titrate the whole mixture with under nitrogen gas bubbling.

※ End point of titration: Appearance of the same color as that of the blank

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20 : reaction time (min)

F : factor of titration solution (50 mM NaOH)

50 : concentration (mM) of titration solution (50 mM NaOH)

0.5 : the volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition (Figure 4).

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