LIPASE [LPM]

mixture of LP(T-01) and MGLP II(T-117)

- ★ Advantages
- ① High Reactivity
 - ② Thermal stability
 - ③ Stability in solution

Preparation and Specification

Appearance : White to pale brownish amorphous powder, lyophilized
 Specific activity : LP ; 1500-3500 U/mg
 MGLP II ; 12-30 U/mg

Applications for Diagnostic Test

The enzyme is useful for enzymatic determination of **triglyceride** when coupled with glycerophosphate oxidase (T-60) and glycerol kinase (T-64).


Fig.1 Thermal stability

Stored in R2 of TG reagent
Residual activity (%) after 30 min. at each temperature
Activity assay method: 37°C


Fig.2-1 Dilution linearity using human serum

Sample : Human serum
Analyzer : Hitachi 7150


Fig.2-2 Dilution linearity using human serum

Sample : Human serum
Analyzer : Hitachi 7150


Fig.3-1 Liquid stability

Residual activity (%) after 45 days at 40°C in R2 of TG reagent



Fig.3-2 Liquid stability(Reactivity to TG)


Reactivity after 3, 4 and 5 weeks at 40°C in R2 of TG reagent. (Control : Reagent stored at 4°C)

Assay

■ Principle

The assay is based on the titration of fatty acids liberated in the following reactions:

■ Unit definition

One unit is defined as the amount of enzyme which liberates 1 μ mole of fatty acid per minute at 37°C under the conditions specified in the assay procedure.

■ Reagents

1. Substrate suspension (Olive oil and Adekatol SO-120) 50 g of Olive oil: (Japanese Pharmacopoeia grade) and 50 g of Adekatol SO-120 are suspended with 150 ml of distilled water.
2. Reaction stopper
Mixture of ethanol and acetone (1:1)
3. Indicator
1% (W/V) Phenolphthalein-ethanol solution
4. Titration solution
50 mM NaOH solution
5. Enzyme dilution buffer
0.1 M KH_2PO_4 -NaOH buffer, pH 8.0 containing 0.1% (W/V) BSA and 0.1% (W/V) NaN_3
6. Reagents
Olive oil: (Japanese Pharmacopoeia grade)
Ethanol: (Japanese Pharmacopoeia grade)
Adekatol SO-120: ADEKA CORPORATION
BSA: Millipore Fraction V pH5.2 #81-053

■ Enzyme solution

Accurately weigh about 10 mg of the sample and add enzyme dilution buffer to make a total of 50 ml. Dilute it with enzyme dilution buffer to adjust the concentration to within 2-4 U/ml.

■ Procedure

1. Pipette accurately 5 ml of substrate suspension and 2 ml of distilled water into a test tube (24 mm i.d. \times 200 mmH) and mix to start the preincubation at 37°C.

2. After 10 min, add 0.5 ml of enzyme solution and mix to start the reaction.
※ In the case of a test blank, add 0.5 ml of enzyme dilution buffer in place of enzyme solution.
3. After 20 min, stop the reaction with 16 ml of reaction stopper.
4. Add 3 drops of indicator and titrate the whole mixture with under nitrogen gas bubbling.
※ End point of titration: Appearance of the same color as that of the blank

$$\begin{array}{l} \text{Titration volume} \quad \text{sample} : V_s \\ \quad \quad \quad \quad \quad \text{blank} : V_c \\ \Delta V = (V_s - V_c) \leq 2.5 \text{ ml} \\ V_c \leq 0.6 \text{ ml} \end{array}$$

■ Calculation

$$\text{Activity (U/mg of powder)} = \frac{\Delta V \times F}{20} \times 50 \times \frac{1}{0.5} \times \frac{1}{X}$$

20 : reaction time (min)

F : factor of titration solution (50 mM NaOH)

50 : concentration (mM) of titration solution (50 mM NaOH)

0.5 : the volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition.

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