# GLUTAMINE SYNTHETASE [GS II]

from *Bacillus* sp.

(L-Glutamate: ammonia Ligase (ADP-forming), EC 6.3.1.2)

## Preparation and Specification

Appearance : White to pale brownish lyophilized powder  
 Specific activity : More than 15U/mg  
 Contaminants : NADH oxidase Less than 0.05% (U/U)

## Properties

Michaelis constant	: NH <sub>4</sub> Cl	0.05 mM	
	Glutamic acid	10.9 mM	
Isoelectric point	: pH 4.85		
Molecular weight	: 500KDa (TSK G3000SW <sub>XL</sub> )		
	53KDa (SDS-PAGE)		
Optimum pH	: 8-9		Figure 1
pH stability	: 5.0-9.5 (37°C, 60min)		Figure 2
Optimum temperature	: 60°C (pH8.0, 50mM Tris-HCl buffer)		Figure 3
Thermal stability	: Stable at 60°C and below (pH8.0, 10min)		Figure 4
Storage stability	: -20°C		
Effect of metal ions	: See Table 1		

## Applications for Diagnostic Test

This enzyme is useful for **elimination of ammonia** in a diagnostic reagent.





### GS II

Table 1. Effect of metal ions GS II activity

Metal	%	Metal	%
None	100	AgCl	82
NaCl	101	CsCl	90
KCl	97	NH <sub>4</sub> Cl	69
EDTA	0	ZnSO <sub>4</sub>	0
LiCl	92	CuSO <sub>4</sub>	0


※ 10mM Metal ion with 5mM MgCl<sub>2</sub>

## Assay

### Principle

The assay is based on the decrease in absorbance at 340 nm as the consumption of NADH proceeds in the following reactions:

PK : Pyruvate kinase  
LDH : Lactate dehydrogenase  
PEP : Phosphoenol pyruvate

### Unit definition

One unit is defined as the amount of enzyme which produces 1  $\mu$ mole of glutamine per minute at 37 °C under the conditions specified in the assay procedure.

### Reagents

1. Reaction mixture		
0.2 M Tris-HCl buffer pH 8.0		0.75 ml
60 mM ATP solution pH 8.0		0.30 ml
1 M NH <sub>4</sub> Cl solution		0.30 ml
30 mM PEP solution		0.20 ml
7.86 mM NADH solution		0.10 ml
0.6 M MgCl <sub>2</sub> solution		0.15 ml
100 U/ml PK <sup>1)</sup> solution		0.10 ml
340 U/ml LDH <sup>2)</sup> solution		0.10 ml

0.45 M Glutamic acid Na	0.60 ml
Distilled water	0.40 ml

1): 100U/ml PK

Stock solution (2000U/ml) of PK with distilled water to make 20-fold solution.

2): 340U/ml LDH

Stock solution (5500U/ml) of LDH with distilled water to make 16-fold solution.

2. Enzyme dilution buffer

10mM Tris-HCl buffer pH8.0

3. Reagent

Tris (hydroxymethyl) aminomethane:

Sigma Chemical Co. #T-1503

ATP (2Na·3H<sub>2</sub>O·Reduced form): Kyowa Hakko Co., Ltd

NADH (2Na·3H<sub>2</sub>O·Reduced form): Kyowa Hakko Co., Ltd

PEP [Phospho(enol)pyruvic acid

tri(cyclohexylammonium) salt] :

Sigma Chemical Co. #P-7252

NH<sub>4</sub>Cl:FUJIFILM Wako Pure Chemical Corporation

Special grade #015-02991

MgCl<sub>2</sub>·6H<sub>2</sub>O:

FUJIFILM Wako Pure Chemical Corporation

Special grade #131-00162

PK:Roche Diagnostics GmbH #10 128 163 001

LDH:Roche Diagnostics GmbH #10 127 876 001

Glutamic Acid Monosodium Salt:

FUJIFILM Wako Pure Chemical Corporation

Special grade #198-02035

### Enzyme solution

Accurately weigh about 10 mg of the sample and add enzyme dilution buffer to make a total of 10 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

## ■ Procedure

1. Pipette accurately 3.0 ml of reaction mixture into a small test tube and Preincubate at 37°C
2. After 5min. add accurately 60  $\mu$ l of enzyme solution and mix to start the reaction at 37°C  
※ In the case of a test blank add 60  $\mu$ l of enzyme dilution buffer in place of enzyme solution
3. After starting the reaction, measure the rate of increase per minutes in absorbance at 340 nm. The rate must be measured within the linear portion of the absorbance curve. (Ex. Linear range from 3 min to 6 min)

Absorbance sample : As/min  
blank : Ab/min

$$\Delta A/\text{min} = (As/\text{min} - Ab/\text{min}) \leq 0.15 \text{ Abs/min}$$

## ■ Calculation

$$\text{Activity (U/mg)} = \frac{\Delta A/\text{min}}{6.3} \times \frac{3.06}{0.06} \times \frac{1}{X}$$

6.3 : millimolar extinction coefficient of NADH at 340nm  
( $\text{cm}^2 / \mu\text{mol}$ )

3.06 : final volume (ml)

0.06 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution

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