(Diagnostic Reagent Grade)


ASAHI KASEI ENZYMES T-39 T-138

PHOSPHOLIPASE D [PLDP (T-39), PLDPV (T-138)]

(Glycerophospholipid specific)

from *Streptomyces* sp.

(Phosphatidyl choline phosphatidohydrolase, EC 3.1.4.4)

Preparation and Specification

PLDP (T-39)

Appearance : Light brownish amorphous powder, lyophilized
Specific activity : More than 100 U/mg solid

PLDPV (T-138)

Appearance : White to brownish lyophilized powder
Specific activity : More than 30 U/mg solid
※ Animal derived material free.


Properties

Substrate specificity	: See Table 1	
Molecular weight	: 46 kDa (gel filtration)	
Isoelectric point	: pH 4.2	
Optimum pH	: 5.5	Figure 1
pH stability	: 4.2-8.5 (37°C, 60 min, 0.05% BSA)	Figure 2
Storage stability	: At least one year at -20°C	Figure 3
Effect of various chemicals	: See Table 2 and Table 3	


Transphosphatidylation Catalyzed by Phospholipase D

Figure 4
Figure 5
Figure 6
Figure 7


Fig.5 Optimum pH of PLDP-Catalyzed

Donor : Uridine 250 mM
 pH6-7 : MES-NaOH buffer
 pH4-6 : Acetate buffer

pH3.5-4 : Glycine-HCL
 pH4-5.5 : Acetate buffer
 Cytidine, Deoxycytidine : 250mM

pH3-4 : Glycine-HCL
 pH4-6 : Acetate buffer
 Adenosine : 50mM
 Deoxyadenosine : 50mM

Phospholipid acceptor : Dipalmitoylphosphatidylcholine in Chloroform (1ml)
 Buffer : 0.5ml contained 1.8 Units of PLDP
 Reaction Temp : 35°C
 Reaction period : 90 min.

Fig.6 The Influence of Reaction Temperature on the Yield of 5'-Phosphatidylcholine by PLDP-Catalyzed Transphosphatidylation

(Donor:dipalmitoylphosphatidylcholine, Acceptor: Cytidine, pH4.5)



Fig.7 The Effect of Acceptor-Donor Ratio on the Yield of 5'-Phosphatidylcytidine by PLDP-Catalyzed Transphosphatidylation




Organic solvent : 10 mM Dipalmitoylphosphatidylcholine in Chloroform (1ml)
 Aqueous solvent : 0.1 M Glycine-HCL (pH4.5) 0.5ml containing 1.8 Units of PLDP
 Reaction : at 35°C, with continuous stirring, for 180min.

Assay

Principle

The assay is based on the increase in absorbance at 500 nm as the formation of quinoneimine dye proceeds in the following reactions:

COD: Choline oxidase

Unit definition

One unit is defined as the amount of enzyme which hydrolyzes 1 μmole of phosphatidylcholine to phosphatidic acid and choline per minute at 37°C under the conditions specified in the assay procedure.

■ Reagents

- Reaction mixture for the first reaction

0.2 M DMG-NaOH buffer pH 5.5	0.10 ml
10 mM CaCl ₂	0.05 ml
3% (W/V) Triton X-100 solution	0.05 ml
10 mM Substrate solution ¹⁾	0.15 ml
Distilled water	0.15 ml

 DMG: 3, 3'-Dimethylglutarate
 1): 10 mM Substrate solution
 Dissolve 78.6 mg of L- *a* -phosphatidylcholine, dioleoyl with 10 ml of 5% Triton X-100 (W/V) .
- Reaction mixture for the second reaction

15 mM 4-AA	0.05 ml
0.2% (W/V) Phenol	0.05 ml
1 M Tris-HCl buffer, pH 8.0	0.05 ml
50 U/ml POD ²⁾	0.05 ml
50 U/ml COD ³⁾	0.05 ml
Distilled water	0.25 ml

 2) : 50 U/ml POD
 Dissolve 500 U (PPU) of POD with 10 ml of distilled water.
- 50 U/ml COD
 Dissolve 500 U of COD with 10 ml of 10 mM Tris-HCl buffer, pH 8.0.
- Reaction stopper
 1 M Tris-HCl buffer containing 10 mM EDTA and 1% (W/V) Cetyltrimethylammonium chloride
 EDTA: Ethylenediamine tetraacetic acid
- Reaction dilution solution
 1% (W/V) Triton X-100
- Enzyme dilution buffer
 10 mM DMG-NaOH buffer, pH 5.5 containing 0.1% Triton X-100
- Reagents:

DMG : Tokyo Kasei Kogyo Co., Ltd. #D1322
 Triton X-100: Yoneyama Yakuhin Kogyo Co., Ltd. #03821
 L- *a* -Phosphatidylcholine, dioleoyl (C18:1, [cis] -9) : Sigma Chemical Co. #P-6354
 EDTA (2Na · 2H₂O) : SEKISUI MEDICAL Co., Ltd. Special grade
 COD: Asahi Kasei Pharma Corporation #T-05
 4-AA: NACALAI TESQUE, INC. Special grade #01907-52
 Cetyltrimethylammonium chloride: Wako Pure Chemical Industries, Ltd. #087-06032
 POD: Sigma Chemical Co. Type II #P-8250

■ Enzyme solution

Weigh about 20 mg of test sample exactly and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

■ Procedure

- Pipette 0.45 ml of reaction mixture for the first reaction into a small test tube and preincubate them at 37°C.
- After 5 min, add exactly 50 μl of enzyme solution and mix to start the first reaction.

※ In the case of a test blank, add 50 μl of enzyme dilution buffer in place of enzyme solution at this point.

- After 10 min, add 0.50 ml of reaction stopper and mix. On stopping the first reaction, add 0.50 ml of the reaction mixture for the second reaction immediately to start the second reaction at 37°C.
- After 20 min, add 1.5 ml of reaction dilution solution and mix.
- After 10 min, measure the absorbance at 500 nm.

Absorbance	sample : As/min
	blank : Ab/min

$$\Delta A/\text{min} = (As/\text{min} - Ab/\text{min}) \leq 0.40 \text{ Abs/min}$$

■ Calculation

$$\text{Activity (U/mg)} = \frac{\Delta A/10}{12.2 \times 1/2} \times \frac{1}{2} \times \frac{3.00}{0.05} \times \frac{1}{X}$$

12.2 : millimolar extinction coefficient of quinoneimine dye at 500 nm (cm² / μmole)

2 : the multiplier derived from the fact that 1 mole of phosphatidyl choline produces 2 mole of H₂O₂

1/2 : the multiplier derived from the fact that 2 mole of H₂O₂ produces 1 mole of quinoneimine dye

10 : reaction time (min)

3.00 : final volume (ml)

0.05 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition (Figure 3).

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