DIAPHORASE (NADH) [DII]

from *Bacillus megaterium*
(NADH: acceptor oxidoreductase, EC 1.6.5.2)

A : Hydrogen acceptor

Preparation and Specification

Appearance : Yellowish amorphous powder, lyophilized
Specific activity : More than 70 U/mg solid

Properties

Substrate specificity	: See Table 1	
Molecular weight	: 50 kDa (SDS-PAGE) , 160 kDa (gel filtration)	
Isoelectric point	: pH 4.2	
Michaelis constant	: NADH $5.5 \times 10^{-4}\text{M}$	
Optimum pH	: 8.0-9.0	Figure 1
pH stability	: 6.0-9.0 (37°C, 60 min)	Figure 2
Optimum temperature	: 40-45°C	Figure 3
Thermal stability	: Stable at 50°C and below (pH 8.0, 10 min)	Figure 4





Applications for Diagnostic Test

This enzyme is useful for enzymatic determination of **reduced NAD**.

Table 1. Substrate specificity


Substrate	Relative activity (%)
NADH	100
NADPH	16

Assay

Principle

The assay is based on the increase in absorbance at 550 nm as the formation of formazane dye (NTBH₂) proceeds in the following reaction:

NADH: Nicotinamide adenine dinucleotide

NTB: Nitrotetrazolium blue

Unit definition

One unit is defined as the amount of enzyme which oxidizes 1 μmole of NADH to NAD⁺ per minute at 37°C under the conditions specified in the assay procedure.

Reagents

- Reaction mixture

0.2 M KH ₂ PO ₄ -NaOH buffer pH 8.0	0.50 ml
0.25% (W/V) NTB solution	0.10 ml
1% (W/V) BSA solution	0.10 ml
10 mM NADH solution	0.10 ml
Distilled water	0.20 ml
- Reaction stopper
 - 0.1 N HCl solution
- Enzyme dilution buffer
 - 0.1 M KH₂PO₄-NaOH buffer pH 8.0 containing
 - 0.1% (W/V) BSA
- Reagents
 - NTB: Dojindo Laboratories # 344-02033
 - BSA: Millipore Fraction V pH5.2 #81-053
 - NADH (2Na·3H₂O·reduced form):

Kyowa Hakko Co. Ltd.

Enzyme solution

Accurately weigh about 20 mg of the sample and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

Procedure

- Pipette accurately 1.0 ml of reaction mixture into a small test tube and preincubate at 37°C.
- After 5 min, add exactly 100 μl of enzyme solution and mix to start the reaction at 37°C.
- At 10 min after starting the reaction, add 2.0 ml of the reaction stopper to stop the reaction.
- Measure the absorbance at 550 nm.

$$\begin{aligned} \text{Absorbance sample} &: A_s \\ \text{Absorbance blank} &: A_b \\ \Delta A &= (A_s - A_b) \leq 0.370 \text{ Abs} \end{aligned}$$

Calculation

$$\text{Activity (U/mg of powder)} = \frac{\Delta A / 10}{12.4} \times \frac{3.10}{0.10} \times \frac{1}{X}$$

12.4 : millimolar extinction coefficient of Formazane dye at 550 nm (cm²/ μmole)

10 : reaction time (min)

3.10 : final volume (ml)

0.10 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition.

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