
Table 1. Substrate specificity


Substrate (1 mM)	Relative activity (%)	Km (mM)
Tyramine	100	0.027
Monoethanolamine	5	15
Trimethylamine	0	
Triethanolamine	0	
Glutamine	0	
Ethanediamine	0	
Dopamine	33	0.23
Histamine	68	0.48
Agmatine	9	0.26
Amylamine	132	0.20
Methylamine	0	
Phenylethylamine	141	0.003
Propylamine	0	
Arylamine	14	0.13
Benzylamine	0	
Hexylamine	115	0.10
1,4-Diaminobutane	0	
Cadaverine	0	

Fig.1 pH Optimum


○ : 3,3-Dimethylglutarate-NaOH buffer  
 △ : Phosphate buffer  
 ● : Tris-HCl buffer  
 ▲ : Glycine-NaOH buffer

Fig.2 pH Stability

37°C, 60 min.  
 ○ : 3,3-Dimethylglutarate-NaOH buffer  
 △ : Phosphate buffer  
 ● : Tris-HCl buffer

Fig.3 Thermal Stability


pH 7.5, 5 min.  
 20 mM Phosphate buffer

## Assay

### Principle

The assay is based on the increase in absorbance at 480 nm as the formation of quinoneimine dye proceeds in the following reactions:

### Unit definition

One unit is defined as the amount of enzyme which oxidizes 1  $\mu$ mole of tyramine to *p*-hydroxybenzylaldehyde per minute at 37°C under the conditions specified in the assay procedure.

## Reagents

### 1. Reaction mixture

0.1 M KH <sub>2</sub> PO <sub>4</sub> -NaOH buffer pH 7.5	0.20 ml
0.2% (W/V) Phenol solution	0.10 ml
15 mM 4-AA solution	0.10 ml
100 U/ml POD solution <sup>1)</sup>	0.05 ml
1 mM Tyramine solution	0.10 ml
Distilled water	0.45 ml

### 1): 100 U/ml POD solution

Dissolve 1,000 U (PPU) with 10 ml of distilled water.

### 2. Reaction stopper

Ethanol (not diluted)

### 3. Enzyme dilution buffer

20 mM KH<sub>2</sub>PO<sub>4</sub> -NaOH buffer pH 7.5

### 4. Reagents

Tyramine: Tokyo Kasei Kogyo Co., Ltd. # A0302

4-AA: NACALAI TESQUE, INC. Special grade

# 01907-52

POD: Sigma Chemical Co. Type II # P-8250

**Enzyme solution**

Accurately weigh about 20 mg of the sample and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

**Procedure**

- Pipette accurately 1.0 ml of reaction mixture into a small test tube and preincubate at 37°C.
- After 5 min, add exactly 50  $\mu$ l of enzyme solution and mix to start the reaction at 37°C.

※ In the case of a test blank, add 50  $\mu$ l of enzyme dilution buffer in place of enzyme solution.

- At 10 min after starting the reaction, add 2.0 ml of the reaction stopper to stop the reaction.
- Measure the absorbance at 480 nm.

Absorbance sample : As  
blank : Ab

$$\Delta A = (A_s - A_b) \leq 0.30 \text{ Abs}$$

**Calculation**

$$\text{Activity (U/mg of powder)} = \frac{\Delta A/10}{17.17 \times 1/2} \times \frac{3.05}{0.05} \times \frac{1}{X}$$

- 17.17 : millimolar extinction coefficient of quinoneimine dye at 480 nm ( $\text{cm}^2 / \mu\text{mole}$ )  
 1/2 : a multiplier derived from the fact that 2 mole of  $\text{H}_2\text{O}_2$  produces 1 mole of quinoneimine dye  
 10 : reaction time (min)  
 3.05 : final volume (ml)  
 0.05 : volume of enzyme solution (ml)  
 X : concentration of the sample in enzyme solution (mg/ml)

**Storage**

Storage at -20°C in the presence of a desiccant is recommended.

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