
(Diagnostic Reagent Grade)

INQUIRY NEEDED
ASAHI KASEI ENZYMES T-33

GLYCEROPHOSPHORYLCHOLINE PHOSPHODIESTERASE [GPCP]

from *Gliocladium roseum*

(sn-Glycero-3-phosphocholine glycerophosphohydrolase, EC 3.1.4.2)

Preparation and Specification

Appearance : White to brownish amorphous powder, lyophilized
 Specific activity : More than 13 U/mg solid

Properties

Substrate specificity : See Table 1
 Isoelectric point : pH 3.75
 Michaelis constant : Glycerophosphorylcholine $1.0 \times 10^{-4}\text{M}$
 Optimum pH : 8.5-9.0
 pH stability : 7.0-9.0 (75°C, 60 min)
 Thermal stability : Stable at 50°C and below (pH 8.0, 10 min)
 Effect of various chemicals : See Table 2 and Table 3
 Activator : Ca^{2+}
 Inhibitors : Zn^{2+} , EDTA

Figure 1
Figure 2
Figure 3

Applications for Diagnostic Test

This enzyme is useful for enzymatic determination of **phospholipid** when coupled with lysophospholipase (T-32) and choline oxidase (T-05).


Table 1. Substrate specificity of GPCP

Substrate	Relative activity (%)
Phosphatidylcholine	0
Lysophosphatidylcholine	0
Glycerophosphorylcholine	100
L- α -Glycerophosphate	0


Table 2. Effect of detergents on GPCP activity

Detergent	Concentration (%)	Relative activity (%)
None	0	100
Triton X-100	0.5	104
Adekator SO-145	0.5	108
Adekator PC-8	0.5	101
Adekator NP-700	0.5	103
Pluronic L-61	0.5	102
Sodium lauryl benzene sulfonate	0.5	0
Sodium lauryl sulfate	0.5	0
Cethyl pyridinium chloride	0.5	77
Cethyl trimethylammonium chloride	0.5	61

Table 3. Effect of metal ions on GPCP activity


Metal ion	Concentration (mM)	Relative activity (%)
None	0	100
CaCl ₂	1	117
MgCl ₂	1	103
MnCl ₂	1	16
BaCl ₂	1	41
CoCl ₂	1	6
CuCl ₂	1	6
NiCl ₂	1	0
NH ₄ Cl	100	115
LiCl	100	112
NaCl	100	109
KCl	100	112
EDTA	1	0

Fig.1 pH Optimum


○ : 3,3-Dimethylglutarate-NaOH buffer
 △ : Phosphate buffer
 ● : Tris-HCl buffer
 ▲ : Glycine-NaOH buffer

Fig.2 pH Stability

37°C, 60 min.
 ○ : 3,3-Dimethylglutarate-NaOH buffer
 ● : Tris-HCl buffer

Fig.3 Thermal Stability

pH 8.0, 10 min.
 10 mM Tris-HCl buffer

Assay


Principle

The assay is based on the increase in absorbance at 500 nm as the formation of quinoneimine dye proceeds in the following reactions:


GPCP

COD

POD

GPCP: Glycerophosphorylcholine

COD: Choline oxidase

Unit definition

One unit is defined as the amount of enzyme which hydrolyzes 1 μmole of glycerophosphorylcholine per minute at 37°C under the conditions specified in the assay procedure.

Reagents

1. Reaction mixture for the first reaction

0.2 M Tris-HCl buffer pH 8.0	0.10 ml
10 mM GPC solution ¹⁾	0.10 ml
10 mM CaCl ₂ solution	0.05 ml
Distilled water	0.25 ml

1): 10 mM GPC solution pH 8.0

Accurately weigh 100 mg of GPC into a brown test tube and add 4 ml of distilled water and 80 μl of 5 N

NaOH. After suspending it, centrifugate at 3,000 rpm for 10 min. Remove and store supernatant. Add 4 ml of distilled water to the precipitate and suspend again. Remove and store supernatant after centrifugation.

Combine the supernatant portions and adjust pH to 8.0 (25°C) with diluted HCl and add distilled water to make a total of 20 ml.

2. Reaction mixture for the second reaction

0.2 M Tris-HCl buffer pH 8.0	0.10 ml
0.1 M EDTA solution pH 8.0	0.20 ml
0.3% 4-AA solution	0.10 ml
0.2 % (W/V) Phenol solution	0.10 ml
60 U/ml COD solution ²⁾	0.10 ml
100 U/ml POD solution ³⁾	0.05 ml
Distilled water	0.35 ml

EDTA: Ethylenediaminetetraacetic acid

2): 60 U/ml COD solution

Dissolve 600 U of COD with 10 ml of 10 mM Tris-HCl buffer pH 8.0.

3): 100 U/ml POD solution

1,000 U (PPU) of POD with 10 ml of distilled water.

3. Enzyme dilution buffer

10 mM Tris-HCl buffer pH 8.0 containing 0.05% (W/V) BSA.

4. Reagents

GPC: Sigma Chemical Co. #G-8005

EDTA (2Na·2H₂O):

KISHIDA CHEMICAL Co., Ltd. #060-29133

4-AA: NACALAI TESQUE, INC. Special grade #01907-52

COD: Asahi Kasei Pharma Corporation #T-05

POD: Sigma Chemical Co. Type II #P-8250

BSA: Millipore Fraction V pH 5.2 #81-053

■ Enzyme solution

Weigh about 20 mg of test sample exactly and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

■ Procedure

1. Pipette 0.45 ml of reaction mixture for the first reaction

into a small test tube and preincubate at 37°C.

2. After 5 min, add exactly 50 μ l of enzyme solution and mix to start the first reaction at 37°C.

※ In the case of a test blank, add 50 μ l of enzyme dilution buffer in place of enzyme solution.

3. After 10 min, add 1.0 ml of reaction mixture for the second reaction and mix to start the second reaction at 37°C.

4. After 20 min, add 1.5 ml of distilled water to stop the reaction.

5. Measure the absorbance at 500 nm.

Absorbance sample : As

blank : Ab

$$\Delta A = A_s - A_b$$

■ Calculation

$$\text{Activity (U/mg)} = \frac{\Delta A/10}{12.0 \times 1/2} \times \frac{1}{2} \times \frac{3.00}{0.05} \times \frac{1}{X}$$

12.0 : millimolar extinction coefficient of quinoneimine dye at 500 nm (cm² / μ mole)

2 : a multiplier derived from the fact that 1 mole of glycerophosphorylcholine produces 2 mole of H₂O

1/2 : a multiplier derived from the fact that 2 mole of H₂O₂ produces 1 mole of quinoneimine dye

10 : reaction time (min)

3.00 : final volume (ml)

0.05 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20 °C in the presence of a desiccant is recommended.

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