


Citrate pump unit for Plasauto Σ TM

Plasauto Σ TM-CP



- Gives extra option of Citrate Anticoagulation
- Specifically designed and engineered for Plasauto Σ
- Safe operation is assured with direct connection to Plasauto Σ by an interlocking system



Citrate pump unit for PlasautoΣ

PlasautoΣTM-CP

Specifications

Citrate pump	Flow rate control range	0.1 to 10mL/min
	Accuracy	Less than ±3% (Over 300mL)
Scale	Rating	10kg
	Accuracy	Less than ±5g of the displayed value of 1,000g
Electric shock protection	Class I equipment, Type BF applied part	
Waterproofing	IPX1	
EMC	Conform to EMC standard EN 60601-1-2:2015 +A1:2021	
Dimensions	177(W) x 281(H) x 164(D) mm	
Weight	2.2kg ±10%	



◆ PlasautoΣ-CP



◆ Connected with PlasautoΣ

AsahiKASEI

Global distributor:
ASAHI KASEI MEDICAL CO., LTD.
1-1-2 Yurakucho, Chiyoda-ku, Tokyo 100-0006 Japan
TEL: +81-3-6699-3771 FAX: +81-3-6699-3773
www.asahikasei-medical.com

Authorized Representative in the European Community:
ASAHI KASEI MEDICAL EUROPE GmbH
Herriotstrasse 1, 60528 Frankfurt am Main Germany
TEL: +49 (0) 69 66 37 15 00 FAX: +49 (0) 69 6 66 51 93

Manufacturer:
MED-TECH INC.
1-103-66, Yoshinodai, Kawagoe-shi, Saitama 350-0833 Japan