



Diagnostics Department
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682
diagnostics @om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41D01-85
March 1, 2009
[5th edition]

Product Name: **DIAPHORASE (NADH) [DI]**
(NADH:acceptor oxidoreductase, EC 1. 6. 5. 2)

Product Number: **T-06**

Origin: *Bacillus megaterium*

Description: Yellowish lyophilized powder

Activity: 30~60 U/mg

Appearance of solution: Clear
(2U/mL in buffer solution)
A400nm ≤ 0.050
A500nm ≤ 0.025
A650nm ≤ 0.010

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant