



SPECIFICATION

ISO: 75F41A09-85
March 1, 2009
[4th edition]

Product Name: **ALKALINE PHOSPHATASE [ALP]**
(Orthophosphoric monoester phosphohydrolase, EC 3. 1. 3. 1)

Product Number: **T-08**

Origin: *Escherichia coli*

Description: White lyophilized powder

Activity: ≥ 40 U/mg

Appearance of solution: Clear
(50U/mL in buffer solution)

| | |
|--------|--------------|
| A400nm | ≤ 0.015 |
| A500nm | ≤ 0.015 |
| A650nm | ≤ 0.015 |

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant