



SPECIFICATION

ISO: 75F41G06-85
March 1, 2009
[5th edition]

Product Name: L- α -GLYCEROPHOSPHATE OXIDASE [GPOM]
(L-Glycero-3-phosphate:oxygen oxidoreductase, EC 1. 1. 3. 21)

Product Number: T-107

Origin: *Streptococcus* sp.

Description: Yellowish lyophilized powder

Activity: ≥ 15.0 U/mg

Appearance of solution: Clear
(50U/mL in buffer solution)
A400nm ≤ 0.150
A500nm ≤ 0.060
A650nm ≤ 0.015

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant