



## SPECIFICATION

ISO: 75F41P04-85  
March 1, 2009  
[ 4<sup>th</sup> edition ]

Product Name: **PHOSPHOLIPASE-C [PLC]**  
(Phosphatidylcholine cholinephosphohydrolase, EC 3. 1. 4. 3)

Product Number: **T-11**

Origin: *Bacillus cereus*

Description: White to brownish lyophilized powder

Activity:  $\geq 30$  U/mg

Appearance of solution: Clear  
(10U/mL in buffer solution)  
A400nm  $\leq 0.100$   
A500nm  $\leq 0.050$   
A650nm  $\leq 0.025$

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.  
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

-----  
M.TOBE  
Plant Director  
Ohito diagnostics plant

H. Nagai

-----  
H.NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant