



SPECIFICATION

ISO: 75F41A03-85
March 1, 2009
[4th edition]

Product Name: **ACYL-COA SYNTHETASE [ACS]**
(Acid:CoA Ligase(AMP), EC 6. 2. 1. 3)

Product Number: **T-16**

Origin: *Pseudomonas fragi*

Description: White lyophilized powder

Activity: 2~8 U/mg

Appearance of solution: Clear
(2U/mL in buffer solution)
A400nm ≤0.050
A500nm ≤0.025
A650nm ≤0.012

Contaminating enzymes:
Catalase ≤1.0 (U/U)
Glucose oxidase ≤0.1 % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant