



SPECIFICATION

ISO: 75F41A02-85
March 1, 2009
[4th edition]

Product Name: **ACYL-COA OXIDASE [ACOD]**
(Acyl-CoA: oxygen oxidoreductase, EC 1. 3. 3. 6)

Product Number: **T-17**

Origin: *Arthrobacter* sp.

Description: Yellowish lyophilized powder

Activity: 20~40 U/mg

Appearance of solution: Clear
(20U/mL in buffer solution)
A400nm ≤ 0.080
A500nm ≤ 0.040
A650nm ≤ 0.020

Contaminating enzymes:
Catalase ≤ 1.00 (U/U)
Glucose oxidase ≤ 0.01 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M.TOBE
Plant Director
Ohito diagnostics plant

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant