



SPECIFICATION

ISO: 75F41B01-85
May 20, 2010
【 3rd edition 】

Product Name: **BILIRUBIN OXIDASE [BOD II]**
(bilirubin:oxygen oxidoreductase, EC 1. 3. 3. 5)

Product Number: **T-188**

Origin: *Bacillus subtilis*

Description: Light blue to blue lyophilized powder

Activity: ≥ 5.0 U/mg

Appearance of solution:
(1mg/mL in buffer solution) Clear
A400nm ≤ 0.050
A500nm ≤ 0.025
A650nm ≤ 0.050

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant