



SPECIFICATION

ISO: 75F41C04-85
March 1, 2009
[5th edition]

Product Name: **CHOLESTEROL ESTERASE [CEN]**
(Sterol-ester acylhydrolase, EC 3. 1. 1. 13)

Product Number: **T-18**

Origin: *Pseudomonas* sp.

Description: A pale brownish lyophilized powder

Activity: 100~250 U/mg

Activity/Protein: ≥ 100 U/mgp

Appearance of solution: Clear
(150U/mL in buffer solution)
A400nm ≤ 0.500
A500nm ≤ 0.250
A650nm ≤ 0.125

Contaminating enzymes:
Catalase ≤ 1.0 % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant