



SPECIFICATION

ISO: 75F41G04-85
March 1, 2009
[2nd edition]

Product Name: **GLYCEROPHOSPHORYLCHOLINE
PHOSPHODIESTERASE [GPCP]**
(sn-Glycero-3-phosphocholine glycerophosphohydrolase, EC 3. 1. 4. 2)

Product Number: **T-33**

Origin: *Gliocladium roseum*

Description: White to brownish lyophilized powder

Activity: ≥ 10.0 U/mg

Appearance of solution: Clear
(50U/mL in buffer solution) A400nm ≤ 0.100
A500nm ≤ 0.050
A650nm ≤ 0.050

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant