



SPECIFICATION

ISO: 75F41A08-85
March 1, 2009
[4th edition]

Product Name: **ALCOHOL OXIDASE [ALOD]**
(Alcohol: oxygen oxidoreductase, EC 1. 1. 3. 13)

Product Number: **T-38**

Origin: *Candida* sp.

Description: Light yellowish lyophilized powder

Activity: 7~20 U/mg

Appearance of solution: Clear
(12U/mL in buffer solution)
A400nm \leq 0.200
A500nm \leq 0.100
A650nm \leq 0.030

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant