



SPECIFICATION

ISO: 75F41H04-85
March 1, 2009
[4th edition]Product Name: **3 α -HYDROXYSTEROID DEHYDROGENASE [3 α -HSD II]**(3 α -Hydroxysteroid:NAD oxidoreductase, EC 1. 1. 1. 50)

Product Number:	T-58
Origin:	<i>Pseudomonas</i> sp.
Description:	White lyophilized powder
Activity:	≥ 30 U/mg
Appearance of solution: (5U/mL in buffer solution)	Clear A400nm ≤ 0.010 A500nm ≤ 0.010 A650nm ≤ 0.010
Contaminating enzymes: NADHoxidase	≥ 90 % (NADH Residual)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M.TOBE
Plant Director
Ohito diagnostics plant-----
H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant