



Diagnostics Department
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682
diagnostics @om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41H02-85
March 1, 2009
[4th edition]

Product Name: **3-HYDROXYBUTYRATE DEHYDROGENASE [3-HBDH II]**

(D-3-Hydroxybutyrate:NAD oxidoreductase, EC 1. 1. 1. 30)

Product Number:	T-68
Origin:	<i>Alcaligenes faecalis</i>
Description:	White lyophilized powder
Activity:	≥1500 U/mg
Appearance of solution: (1mg/mL in buffer solution)	Clear A400nm ≤0.050 A500nm ≤0.040 A650nm ≤0.030
Contaminating enzymes: NADHoxidase	≤0.0004 % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant