

SPECIFICATION

ISO: 75F41C13-85
March 1, 2009
[4th edition]

Product Name: **CHOLESTEROL OXIDASE [CON II]**
(Cholesterol:oxygen oxidoreductase, EC 1. 1. 3. 6)

Product Number: **T-84**

Origin: *Microorganism*


Description: Yellowish solution

Activity: ≥ 500 U/mL


Appearance of solution: Clear
(5U/mL in buffer solution)
A400nm ≤ 0.030
A500nm ≤ 0.020
A650nm ≤ 0.010

Contaminating enzymes:
Catalase ≤ 2.0 % (U/U)
Cholesterol esterase ≤ 0.02 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -80°C .
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



M.TOBE
Plant Director
Ohito diagnostics plant



H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant