



SPECIFICATION

ISO: 75F41A05-85
March 1, 2009
[5th edition]

Product Name:	ADP-HEXOKINASE [ADP-HKP II] (D-Hexose-6-phosphotransferase, EC 2. 7. 1. 147)
Product Number:	T-92
Origin:	<i>Pyrococcus furiosus</i>
Description:	White lyophilized powder
Activity:	≥ 30 U/mg
Appearance of solution: (1mg/mL in buffer solution)	Clear A400nm ≤ 0.015 A500nm ≤ 0.015 A650nm ≤ 0.015
Contaminating enzymes:	
ATPase	≤ 0.01 % (U/U)
NADPHoxidase	≤ 0.01 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M.TOBE
Plant Director
Ohito diagnostics plant

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant