



SPECIFICATION

ISO: 75F41A06-85
March 1, 2009
[4th edition]

Product Name: **ADP-HEXOKINASE [ADP-HKT II]**
 (D-Hexose-6-phosphotransferase, EC 2. 7. 1. 147)

Product Number: **T-93**

Origin: *Thermococcus litoralis*

Description: White lyophilized powder

Activity: ≥ 25 U/mg

Appearance of solution: Clear
(1mg/mL in buffer solution) A400nm ≤ 0.015
 A500nm ≤ 0.015
 A650nm ≤ 0.015

Contaminating enzymes: ATPase ≤ 0.01 % (U/U)
 NADPHoxidase ≤ 0.01 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant