



Diagnostics Department  
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan  
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682  
diagnostics @om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41A07-85

March 1, 2009

[ 5<sup>th</sup> edition ]

Product Name: 1,5-ANHYDRO GLUCITOL-6-PHOSPHATE DEHYDROGENASE  
[AG6PDH II]

( , EC . . . )

Product Number:	T-95
Origin:	<i>Microorganism</i>
Description:	White lyophilized powder
Activity:	≥20 U/mg
Appearance of solution: (1mg/mL in buffer solution)	Clear A400nm ≤0.015 A500nm ≤0.015 A650nm ≤0.015
Contaminating enzymes:	
ATPase	≤0.05 % (U/U)
NADPHoxidase	≤0.02 % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.  
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

*M. Toke*  
-----  
M.TOBE  
Plant Director  
Ohito diagnostics plant

*H. Nagai*  
-----  
H.NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant