



SPECIFICATION

ISO: 75F41C03-85

March 1, 2009

[4th edition]

Product Name: **CHOLESTEROL ESTERASE [CEBP-M]**
(Sterol-ester acylhydrolase, EC 3. 1. 1. 13)

Product Number: **T-98**

Origin: *Microorganism*

Description: White to off white lyophilized powder

Activity: ≥ 10.0 U/mg

Appearance of solution: Clear
(1mg/mL in buffer solution)
A400nm ≤ 0.030
A500nm ≤ 0.015
A650nm ≤ 0.010

Contaminating enzymes:
Catalase ≤ 1.0 % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant